

AVOIRDUPOIS CAST IRON SLOTTED WEIGHTS

Weight	Part #	Width (diameter) base (inches)	Height base (inches)	Slot Width (inches)	Stem Diameter (inches)	Stem Height (inches)	Hanger Capacity
200 lb*	9752	17.00	3.63	1.00			
100 lb*	9750	17.00	1.88	1.00			
50 lb	9498	10.25	2.75	0.56			
25 lb	9497	10.25	1.38	0.56			
20 lb	9496	8.00	1.88	0.56			
10 lb	9494	6.00	1.50	0.56			
5 lb	9492	5.00	1.13	0.56			
2 lb	9490	5.00	0.56	0.56			
1 lb	9488	5.00	0.31	0.56			
1/2 lb	9486	3.75	0.31	0.56			
5 lb Hanger	9484	5.00	1.00		0.31	36.00	500 lbs
2 lb Hanger	9482	3.75	1.13		0.25	20.25	300 lbs
2 lb Hanger	9480	3.75	1.13		0.25	9.25	300 lbs
1 lb Hanger	9478	3.75	0.31		0.25	20.25	200 lbs
1 lb Hanger	9476	3.75	0.31		0.25	9.25	200 lbs

*Signifies weights with handles