

# OIML WEIGHT SETS CONFIGURATION TABLE

		1 mg	2 mg	5 mg	10 mg	20 mg	50 mg	100 mg	200 mg	500 mg	1 g	2 g	5 g	10 g	20 g	50 g	100 g	200 g	500 g	1 kg	2 kg	5 kg	10 kg	20 kg	
Weight Set	Part #																								
500 mg—1 mg	7440	1	2	1	1	2	1	1	2	1															
5 g—1 mg	7441	1	2	1	1	2	1	1	2	1	1	2	1												
50 g—1 mg	7433	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1									
50 g—1 g	7429										1	2	1	1	2	1									
100 g—1 mg	7428	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1								
100 g—1 g	7424										1	2	1	1	2	1	1								
200 g—1 mg**	7423	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2							
200 g—1 g**	7422										1	2	1	1	2	1	1	2							
200 g—1 mg	7421	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	1							
200 g—1 g	7420										1	2	1	1	2	1	1	1							
500 g—1 mg	7419	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1						
500 g—1 g	7418										1	2	1	1	2	1	1	2	1						
1 kg—1 mg	7417	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1					
1 kg—1 g	7416										1	2	1	1	2	1	1	2	1	1					
2 kg—1 mg*	7415	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2				
2 kg—1 g*	7414										1	2	1	1	2	1	1	2	1	1	2				
2 kg—1 mg	7413	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	1				
2 kg—1 g	7412										1	2	1	1	2	1	1	2	1	1	1				
5 kg—1 mg	7411	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1			
5 kg—1 g	7410										1	2	1	1	2	1	1	2	1	1	2	1			
5 kg—1 kg	7442																			1	2	1			
10 kg—1 mg	7409	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1		
10 kg—1 g	7408										1	2	1	1	2	1	1	2	1	1	2	1	1		
20 kg—1 mg	7407	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	2	1	1	1	
20 kg—1 g	7406										1	2	1	1	2	1	1	2	1	1	2	1	1	1	

\*Weight Set contains two 2 kg weights

\*\*Weight Set contains two 200 g weights