

PIPETTE CALIBRATION SERVICES

	Level 1	Level 2	Level 3
NVLAP accredited	•	•	•
turnaround time	72 hours	72 hours	72 hours
"As found" data	10-10-10 using 3 volumes for each channel	4-4-4 using 3 volumes for each channel	4-4-4 using 3 volumes for each channel
"As left" data	10-10-10 using 3 volumes for each channel	4-4-4 using 3 volumes for each channel	N/A
Thorough pipette cleaning	•	•	
Adjustments as required	•	•	
Seals replaced as needed (other parts extra)	•	•	