

PIPETTE CALIBRATION SERVICES

PIPETTE CALIBRATION SERVICES							
SERVICE LEVEL	NVLAP ACCREDITED	TURNAROUND TIME	“AS FOUND” DATA	“AS LEFT” DATA	THOROUGH CLEANING	ADJUSTMENTS AS REQUIRED	
Level 1 Single Channel	●	1 week	10 - 10 - 10	10 - 10 - 10	●	●	
Level 1 Multi Channel	●	1 week	10 - 10 - 10	10 - 10 - 10	●	●	
Level 1A Single Channel	●	1 week	4 - 4 - 4	10 - 10 - 10	●	●	
Level 1A Multi Channel	●	1 week	4 - 4 - 4	10 - 10 - 10	●	●	
Level 2 Single Channel	●	1 week	4 - 4 - 4	4 - 4 - 4	●	●	
Level 2 Multi Channel	●	1 week	4 - 4 - 4	4 - 4 - 4	●	●	
Level 2A Single Channel	●	1 week	4 - 4	4 - 4	●	●	
Level 2A Multi Channel	●	1 week	4 - 4	4 - 4	●	●	
Level 3 Single Channel	●	1 week	4 - 4 - 4	N/A			
Level 3 Multi Channel	●	1 week	4 - 4 - 4	N/A			