

TROEMNER PIPETTES UNCERTAINTY TABLE

NVLAP CODE 105013-0

Test Volume in μl	Gravimetric Method Best Uncertainty (\pm) in μl	Spectrophotometric method Best Uncertainty (\pm) in μl
0.1	0.027	0.031
0.2	0.025	0.026
0.5	0.025	0.050
1.0	0.025	0.078
2.5	0.020	0.119
5.0	0.089	0.251
10.0	0.128	0.339
25.0	0.14	0.39
50.0	0.30	1.87
100.0	0.34	2.60
500.0	0.78	13.12
1000.0	1.38	16.76
2500.0	5.25	23.08
5000.0	5.98	-