



Troemner Custom Designed Carts to Transport and Protect your Weights!

Troemner manufactures custom carts to meet your needs. This is one example of a custom cart that we did for a customer who needed a cart to transport and protect its weights in a harsh environment. We decided, in this case, to use the 304 stainless steel with a welded stainless steel construction for its strength and protection. The stainless steel lid was designed with two handles and is hinged to make it easy to open and close. The wrap over design of the lid prevents water from entering, which is a very important feature in a wash down, clean room, corrosive or harsh environment. The lid when closed is secured by two stainless steel latches. Inside this cart, each 20 kg or 25 kg weight has its own lined pocket to prevent movement and possible scratching of the weight. Troemner also installed 5 inch non-marking wheels. Two of the wheels in this design swivel and lock for movement and safety.



Troemner works with you to design a cart that matches your requirements. Give us a call at 1.800.352.7718 and ask for Jack Schmidt at ext 216 or ask for Customer Service.

Media contact: Brooke Trexler, Marketing Manager, 856.686.4202, btrexler@troemner.com