



Hummingbird Press Release



*Precision & Balance
Matched Only
By Nature.*

New! Troemner Stainless Steel Heavy Capacity Weights, Resistant to Corrosion, Ideal for Clean Rooms and Wash Down Environments!

Troemner's new Stainless Steel Heavy Capacity Weights are an ideal solution for the food, beverage, pharmaceutical, nuclear and fine chemical industries. Manufactured from stainless steel, these weights are highly stable and resistant to corrosion - ideal for clean rooms, wash down environments and as a replacement for your cast iron weights. Troemner's Stainless Steel Heavy Capacity Weights can be supplied with Troemner's NIST Traceable Certificate up to 2000 kg (5000 lbs) or with a NVLAP Weight Calibration Certificate for weights up to 1000 kg (3000 lbs). The tolerances of these weights are adjusted to ANSI/ASTM Class 4, 5, 6 or NIST Class F.



The Troemner Heavy Capacity Weights are manufactured from a stainless steel alloy which possesses closely controlled density, low magnetic susceptibility, good stability and resistance to corrosion and handling damage. The weights range from 20 kg / 75 lb and up, are cylindrical in shape, have a flat even top surface and have a secure eyebolt on the top center of the weight to allow them be moved easily with various types of industrial equipment. All of Troemner Heavy Capacity Weights have a two-piece construction that allows the weight to be adjusted. Custom Stainless Steel Heavy Capacity Weights of any nominal value for 20 kg / 75 lb up to 2000 kg / 5000 lb can be manufactured to meet your specific requirements. Troemner also manufactures a variety of custom containers and carts to properly store and transport your weights to keep them clean and protected.

For more information visit our website at www.troemner.com or contact us at 1.800.352.7718.



201 Wolf Drive, PO Box 87,
Thorofare, NJ 08086-0087 USA
Phone: 856-686-1600
Fax: 856-686-1601
<http://www.troemner.com>
e-mail: troemner@troemner.com