

# AVOIRDUPOIS STAINLESS STEEL SLOTTED WEIGHTS \*

Weight	Part #	Description	Width (diameter) inches	Height inches	Slot
50 lb	12600-S	Slotted Interlocking	7.88	3.88	0.50
25 lb	12602-S	Slotted Interlocking	6.98	2.55	0.50
20 lb	12604-S	Slotted Interlocking	6.88	2.12	0.50
10 lb	12606-S	Slotted Interlocking	4.95	2.07	0.50
5 lb	12608-S	Slotted Interlocking	4.00	1.65	0.50
3 lb	12609-S	Slotted Interlocking	4.00	1.03	0.50
2 lb	12610-S	Slotted Interlocking	4.00	0.72	0.50
1 lb	12612-S	Slotted Interlocking	2.00	1.42	0.38
0.5 lb	12614-S	Slotted Interlocking	1.63	1.17	0.38
0.2 lb	12616-S	Slotted Interlocking	2.00	0.38	0.38
0.1 lb	12618-S	Slotted Interlocking	1.63	0.31	0.38
8 oz	12620-S	Slotted Interlocking	1.63	1.17	0.38
4 oz	12622-S	Slotted Interlocking	1.25	0.64	0.38
2 oz	12624-S	Slotted Interlocking	1.63	0.37	0.38
50 lb	12700-S	Slotted Flat	7.88	3.78	0.50
25 lb	12702-S	Slotted Flat	6.88	2.52	0.50
20 lb	12704-S	Slotted Flat	6.88	2.52	0.50
10 lb	12706-S	Slotted Flat	4.94	1.99	0.50
5 lb	12708-S	Slotted Flat	4.00	1.55	0.50
3 lb	12709-S	Slotted Flat	4.00	0.93	0.50
2 lb	12710-S	Slotted Flat	4.00	0.62	0.50
1 lb	12712-S	Slotted Flat	2.00	1.33	0.38
0.5 lb	12714-S	Slotted Flat	1.63	1.08	0.38
0.2 lb	12716-S	Slotted Flat	1.63	0.45	0.38
0.1 lb	12718-S	Slotted Flat	1.63	0.22	0.38
8 oz	12720-S	Slotted Flat	1.63	1.08	0.38
4 oz	12722-S	Slotted Flat	1.63	0.54	0.38
2 oz	12724-S	Slotted Flat	1.63	0.27	0.38
1 oz	12726-S	Slotted Flat	1.25	0.26	0.38
1/2 oz	12728-S	Slotted Flat	1.25	0.13	0.38
1/4 oz	12730-S	Slotted Flat	1.25	0.06	0.38
1/8 oz	12732-S	Slotted Flat	1.25	0.03	0.38
1/16 oz	12734-S	Slotted Flat	1.25	0.01	0.38
1/32 oz	12736-S	Slotted Flat	1.25	0.01	0.38

**\*When ordering Stainless Steel Slotted Weights, specify the part number and tolerance class desired. Stainless Steel Slotted Weights are available in tolerance classes ANSI/ASTM Class 1 through Class 6.**