

METRIC STAINLESS STEEL SLOTTED WEIGHTS *

Weight	Part #	Description	Width (diameter) inches	Height inches	Slot
25 kg	12800-S	Slotted Interlocking	7.88	4.27	0.50
20 kg	12802-S	Slotted Interlocking	7.88	3.47	0.50
10 kg	12804-S	Slotted Interlocking	6.88	2.31	0.50
5 kg	12806-S	Slotted Interlocking	4.94	2.32	0.50
3 kg	12808-S	Slotted Interlocking	4.94	1.37	0.50
2 kg	12810-S	Slotted Interlocking	4.00	1.54	0.50
1 kg	12812-S	Slotted Interlocking	4.00	0.82	0.50
1/2 kg	12814-S	Slotted Interlocking	4.00	0.46	0.50
500 g	12816-S	Slotted Interlocking	4.00	0.46	0.50
300 g	12817-S	Slotted Interlocking	3.00	0.47	0.38
200 g	12818-S	Slotted Interlocking	2.00	0.71	0.38
100 g	12820-S	Slotted Interlocking	2.00	0.41	0.38
50 g	12822-S	Slotted Interlocking	1.63	0.34	0.38
25 kg	12900-S	Slotted Flat	7.88	4.15	0.50
20 kg	12902-S	Slotted Flat	7.88	3.33	0.50
10 kg	12904-S	Slotted Flat	6.88	2.22	0.50
5 kg	12906-S	Slotted Flat	4.94	2.10	0.50
3 kg	12908-S	Slotted Flat	4.94	1.28	0.50
2 kg	12910-S	Slotted Flat	4.00	1.44	0.50
1 kg	12912-S	Slotted Flat	4.00	0.72	0.50
1/2 kg	12914-S	Slotted Flat	2.00	1.46	0.38
500 g	12916-S	Slotted Flat	2.00	1.46	0.38
300 g	12917-S	Slotted Flat	2.00	0.87	0.38
200 g	12918-S	Slotted Flat	1.63	0.94	0.38
100 g	12920-S	Slotted Flat	1.63	0.47	0.38
50 g	12922-S	Slotted Flat	1.63	0.24	0.38
20 g	12924-S	Slotted Flat	1.25	0.19	0.38
10 g	12926-S	Slotted Flat	1.25	0.09	0.38
5 g	12928-S	Slotted Flat	1.31	0.10	0.25
3 g	12929-S	Slotted Flat	1.30	0.06	0.25
2 g	12930-S	Slotted Flat	1.30	0.04	0.25
1 g	12932-S	Slotted Flat	1.32	0.02	0.25

***When ordering Stainless Steel Slotted Weights, specify the part number and tolerance class desired. Stainless Steel Slotted Weights are available in tolerance classes ANSI/ASTM Class 1 through Class 6.**