

FLOW METER CALIBRATION SERVICES

	Level 1	Level 2	Level 3
Mass Flow Meters & Controllers	•	•	•
Variable volume meters (rotameters)			• ⁺
ISO 17025 compliant	•	•	•
NIST Traceable	•	•	•
Uncertainty data	•	•	•
Turnaround time	Less than 2 weeks	Less than 2 weeks	Less than 2 weeks
Expedited service available	•	•	•
Helium leak and high pressure checks	•	•	Helium leak only
"As found" data	7 points *	7 points *	7 points typically **
Replace seals	•		
Repair parts/labor	Most repair parts & labor	Repair parts not included	Repair parts not included
Resize flow rate	•	•	
"As left" data	7 points *	7 points *	

* 7 point "As Found" / "As Left Data" @ 0%, 25%, 50%, 75%, 100%, 50%, and 0% full scale flow rate.

Troemner can tailor your calibration with additional or custom points for an additional charge.

+ Rotameters number of points dependent on scale