

Humidity Calibration Service Request for Quote/Order Form

(Please use one form per Humidity Device)

Organization: _____ Contact: _____ Phone: _____

FAX TO: 856-686-1601 Sales Department ATTN: _____

For calibration service, send items and include this form to:
Henry Troemner LLC, 201 Wolf Drive, Thorofare, NJ 08086

****Please complete the Troemner Request for Quote/Order Form Cover Page to fax with this form****

Order or **Request for Quote** (Check one)

Level 1 Humidity Calibration Service - Troemner's most comprehensive offering designed to meet stringent regulatory requirements. Service includes an accredited NIST traceable calibration - 3 points "as found", 3 points "as left" (25%, 50% and 75% RH) at three temperatures (20°C, 23°C and 25°C) - includes stated uncertainty and adjustments.

Level 2 Humidity Calibration Service - For laboratories that require an accredited NIST traceable calibration. Service includes: 3 points "as found" and 3 points "as left" (25%, 50% and 75% RH) at a single calibrated temperature (23°C) - includes stated uncertainty and adjustments.

Level 3 Humidity Calibration Service - A basic single point accredited NIST traceable humidity measurement. Service includes 1 point "as found" (25%, 50% or 75% RH) at a single calibrated temperature (23°C) - includes stated uncertainty. Note - no adjustments or "as left" data provided with this service.

Custom Humidity Certification - Troemner can tailor a humidity calibration to meet your needs. Please note your custom parameters below and Troemner will provide a quote.

Relative Humidity Calibration Point(s) % RH: _____ at °C: _____

Dew Point

Humidity Measuring Equipment Device:

Humidity Probe with System Read Out

Humidity Probe Only

Readout Manufacturer: _____ Model #: _____ Serial #: _____

Probe Manufacturer: _____ Model #: _____ Serial #: _____

Recalibration Interval: _____

Shipping Instructions: _____

3-117-BF (Rev 01/05)

For Troemner Use Only

Quoted Price: _____

Quote #: _____

Date: _____

Name: _____

