



## ISO/IEC 17025 Calibration Certificate

201 Wolf Drive • P.O. Box 87 • Thorofare, NJ 08086-0087 • Phone: 856-686-1600 • Fax: 856-686-1601 • www.troemner.com • e-mail: troemner@troemner.com

Page 2 of 2 Pages

### Dew Point / Frost Point Calibration

#### NAME AND ADDRESS OF CUSTOMER

Troemner, LLC  
201 Wolf Dr.  
Thorofare, NJ 08086

Certificate Number: 123456-DewPoint  
Serial Number: 12345678910  
ID #: 12345678910

#### SECTION 8: CALIBRATION DATA

##### AS FOUND / AS LEFT DATA

Calibration Point °C DP / FP	Standard Indicated °C DP / FP	DUT Indicated °C DP / FP	DUT Error °C DP / FP	Tolerance °C DP / FP	Uncertainty °C DP / FP
-18.7	-18.21	-18.30	-0.09	1.50	0.37
18	18.68	18.70	0.02	1.50	0.37
39	39.35	39.40	0.05	1.50	0.37

#### SECTION 9: ADDENDUM

Customer Supplied Tolerance