



## ISO/IEC 17025 Calibration Certificate

201 Wolf Drive • P.O. Box 87 • Thorofare, NJ 08086-0087 • Phone: 856-686-1600 • Fax: 856-686-1601 • www.troemner.com • e-mail: troemner@troemner.com

Page 2 of 2 Pages

### Temperature-System Calibration

#### NAME AND ADDRESS OF CUSTOMER

Troemner, LLC  
201 Wolf Dr.  
Thorofare, NJ 08086

Certificate Number: 123456-Temperature  
Serial Number: 12345678910  
ID #: 12345678910

#### SECTION 8: CALIBRATION DATA

##### AS FOUND DATA

Calibration Point °C	Standard Indicated °C	DUT Indicated °C	DUT Error °C	Tolerance °C	Uncertainty °C
Channel A					
0	0.050	0.030	-0.020	0.100	0.020
40	40.007	39.983	-0.024	0.100	0.020
100	100.134	99.998	*-0.136	0.100	0.020

\*Denotes Out of Tolerance condition

##### As Found Coefficients

R0 9.990412E+01  
Ap -1.898769E-02  
Bp -1.937566E-05  
An -1.891432E-02  
Bn -3.458489E-04

##### AS LEFT DATA

Calibration Point °C	Standard Indicated °C	DUT Indicated °C	DUT Error °C	Tolerance °C	Uncertainty °C
Channel A					
0	0.030	0.033	0.003	0.100	0.020
40	40.005	40.013	0.008	0.100	0.020
100	100.032	100.042	0.010	0.100	0.020

##### Coefficients as Programmed:

R0 9.989358E+01  
Ap -1.914325E-02

#### SECTION 9: ADDENDUM

Customer Supplied Tolerance