




Regulatory Standards Information

UL = cUL = CSA = TUV ...

- UL, CSA and TUV are equivalent regulatory agencies who test to a variety of safety, electrical, and operational standards
- Choice of agency is preference - how easy it is to work with them, costs associated, etc.
- All agencies provide thorough testing and conduct surprise inspection of facilities
- It is less about the agency and more about the standards they test to
- The Standards & Regulations page of our instruction manuals shows that all safety and regulatory standards are met by these alternate agencies.
- CE is a self declared marking that is required to sell into the European Market, and highly preferred in others
- Part of regulatory testing is ensuring that the equipment meets all CE requirements
- **ALL Troemner equipment** is regulatory tested and carries the CE mark

If you require additional information on Troemner's regulatory requirements or any of our products please call +1.856.686.1600 or email troemner_web@troemner.com



Troemner, LLC